

Fits Virtually All Power Chairs

Harmar's AL500 and AL500 Heavy Duty accommodates virtually any power chair. Available in traditional 350 lbs. or heavy-duty 400 lbs. models.

AL500 is sold separately or bundled with the Golden GP160 Power Chair.

Features & Benefits

- Easy, one-switch operation, license plate light, key switch and courtesy light
- (4) self-tensioning Q'Straint retractors; no chair modifications required
- Manual crank backup and license plate mount included
- Swing-Away option available



Optional AL105 Swing-Away



AL500 Universal Power Chair Lift is also available as a bundle with the AL105 swing away, AL500 auto lift and AL050C cover with the Golden GP160 Power Chair

Specifications	AL500	AL500 HD Heavy Duty
Accommodates	Power wheelchairs	Power wheelchairs
Platform	27.5" x 39" aluminum deck	27.5" x 39" aluminum deck
Lifting Capacity	350 lbs. (159 kg)	400 lbs (181 kg)
Installed Weight	75 lbs. (34 kg)	77 lbs. (35 kg)
Securement	(4) Self-tensioning Q'Straint retractors	(4) Self-tensioning Q'Straint retractors
Hitch	AL123 class III hitch adapter and AL107 hitch post; class II optional	AL123HD class III hitch adapter and AL107 hitch post included
Options	Battery pack, Swing-Away, Power wheelchair covers	Battery pack, Swing-Away, Power wheelchair covers
Warranty	3-year transferable	3-year transferable