



The Bruno Elite Indoor Curved Rail Stair lift offers custom precision and quality. Experience the comfort of a smooth ride with soft starts and stops. The stylized design can be customized with multiple upholstery options, and the painted rail with covered gear rack ensure the Elite will complement any home's décor.

# BRUNO ELITE INDOOR CURVED RAIL STAIR LIFTS



## Features

- 400 lbs. weight capacity
- Maximum speed 25 ft. per minute
- Adjustable width armrests and footrest height
- Multiple position seat height
- Armrest rocker switch operation
- Flip up armrests, seat and footrest
- Two wireless call/send remotes
- Offset swivel seat for safe entry/exit
- Track fits within 7" - 8" of wall
- Safety obstruction sensors
- Seat belt
- Made in the USA
- Bruno Stair Lift Gold warranty: five year major components; two-year parts

## Available Options

- Power swivel seat
- Power folding footrest
- Top or bottom park position, extending rail away from the staircase
- Mid-park and charging station for staircases with landings
- Six different upholstery choices
- Larger seat pad: 21.5" wide x 16" deep
- Larger footrest: 12" long x 11" wide



1 Year Limited Warranty  
on New Equipment Installs

An independently owned and operated franchise business.

✓ Stair Lifts

✓ Vertical Platform Lifts

✓ Pool Lifts

✓ Patient Lifts

✓ Ramps

✓ Auto Lifts

✓ Bath Modifications

✓ Commercial

✓ Residential

✓ And More