

Options

- **Straight Rail Lift System**

- **Overhead lift motor** – The Straight Rail Lift System is compatible with both the Liko Multirall and Liko Likorall overhead lift motors. Find more information about the Liko Multirall here [**<hyperlink to product page>**](#) and the Liko Likorall here [**<hyperlink to product page>**](#)
- **System configurations** – Eight different system configurations for optimum customization:

Liko FreeSpan Straight Rail Lift System Configurations							
Weight capacity*	Rail length	Unit height for transport	Unit height for lifting/transfer	Unit width	Unit depth	Unit weight	
551 lbs (250 kg)	59.1" (1500 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	78.3" - 96.1" (1990 mm - 2440 mm)	48.8" (1240 mm)	91.5 lbs (41.5 kg)	
551 lbs (250 kg)	78.7" (2000 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	98" - 115.7" (2490 mm - 2940 mm)	48.8" (1240 mm)	97.7 lbs (44.3 kg)	
551 lbs (250 kg)	98.4" (2500 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	117.7" - 135.4" (2990 mm - 3440 mm)	48.8" (1240 mm)	103.8 lbs (47.1 kg)	
551 lbs (250 kg)	118.1" (3000 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	137.4" - 155.1" (3490 mm - 3940 mm)	48.8" (1240 mm)	110.2 lbs (50 kg)	
551 lbs (250 kg)	137.8" (3500 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	157.1" - 174.8" (3990 mm - 4440 mm)	48.8" (1240 mm)	116.4 lbs (52.8 kg)	
507 lbs (230 kg)	157.5" (4000 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	176.8" - 194.5" (4490 mm - 4940 mm)	48.8" (1240 mm)	122.6 lbs (55.6 kg)	
441 lbs (200 kg)	177.2" (4500 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	196.5" - 214.2" (4990 mm - 5440 mm)	48.8" (1240 mm)	128.7 lbs (58.4 kg)	
441 lbs (200 kg)	196.9" (5000 mm)	66.9" (1700 mm)	88.6" - 100.4" (2250 mm - 2550 mm)	216.1" - 233.9" (5490 mm - 5940 mm)	48.8" (1240 mm)	135.1 lbs (61.3 kg)	

- **Traverse Lift System**

- **Overhead lift motor** – The Traverse Lift System is compatible with both the Liko Multirall and Liko Likorall overhead lift motors. Find more information about the Liko Multirall here [**<hyperlink to product page>**](#) and the Liko Likorall here [**<hyperlink to product page>**](#)

- **UltraTwin Lift System – Straight Rail**

- **System configurations** – Eight different system configurations for optimum customization:

Liko FreeSpan UltraTwin Straight Rail Lift System Configurations							
Weight capacity*	Rail length	Unit height for transport	Unit height for lifting/transfer	Unit width	Unit depth	Width between rails	Unit weight
1102 lbs (500 kg)	59.1" (1500 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	78.3" (1990 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	173 lbs (78.4 kg)
1102 lbs (500 kg)	78.7" (2000 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	98" (2490 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	185 lbs (84 kg)
1102 lbs (500 kg)	98.4" (2500 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	117.7" (2990 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	198 lbs (89.6 kg)
1102 lbs (500 kg)	118.1" (3000 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	137.4" (3490 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	210 lbs (95.4 kg)
1102 lbs (500 kg)	137.8" (3500 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	157.1" (3990 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	223 lbs (101 kg)
1014 lbs (460 kg)	157.5" (4000 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	176.8" (4490 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	235 lbs (106.6 kg)
882 lbs (400 kg)	177.2" (4500 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	196.5" (4990 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	247 lbs (112.2 kg)
882 lbs (400 kg)	196.9" (5000 mm)	66.9" (1700 mm)	89.4" - 101.2" (2270 mm - 2570 mm)	216.1" (5490 mm)	56.7" (1440 mm)	11.8" - 31.5" (300 mm - 800 mm)	260 lbs (118 kg)

*Weight capacity determined by weight capacity of accessory in use