



BRUNO ELITE CURVED STAIRLIFT



The Bruno Elite Indoor Curved Rail Stairlift offers custom precision and quality. Experience the comfort of a smooth ride with soft starts and stops. The stylized design can be customized with multiple upholstery options, and the painted rail with a covered gear rack ensures the Elite will complement any home's decor.





FEATURES

- 400 lbs. weight capacity
- Maximum speed 25 ft. per minute
- Adjustable width armrests and footrest height
- Multiple position seat height
- Armrest rocker switch operation
- Flip-up armrests, seat, and footrest
- Two wireless call/send remotes
- Offset swivel seat for safe entry/exit
- Track fits within 7"- 8" of the wall
- Safety obstruction sensors
- Seat belt
- Made in the USA
- Bruno Stair Lift Gold warranty: five-year major components; two-year parts

AVAILABLE OPTIONS

- Power swivel seat
- Power folding footrest
- Top or bottom park position, extending rail away from the staircase
- Mid-park and charging station for staircases with landings
- Six different upholstery choices
- Larger seat pad: 21.5" wide x 16" deep
- Larger footrest: 12" long x 11" wide

101 Mobility is a leader in providing a full range of accessibility solutions in your local market. Our local, professional stair lift specialists will meet with you at your home to evaluate the appropriate stair lift solution and provide you with a custom quote, including installation.