

Legacy IITM

STAIRLIFT

Staying Home
ELEVATE YOUR LIFE

“ I love my Legacy Stairlift for the freedom & peace-of-mind it gives me! ”



LIFETIME PARTS WARRANTY

Batteries not included.
Lifetime to the original purchaser.



Legacy II[™]

STAIRLIFT



Custom
colored track
and metal
components



STANDARD FEATURES

- **Lifetime Parts Warranty**- Batteries not included
- Folds to within 12.5" of the wall
- Attaches to steps, not the wall
- Narrow, low profile track (2.5" tall x 6.5" wide)
- Footrest obstruction sensor with audio beep notification
- Battery power allows up to 20 trips during power outages
- 3' Smart technology charge strip at top & bottom of track
- Rocker switch conveniently located on either armrest
- Constant pressure controls
- Wireless remote controls for upper & lower landings
- Chassis mounted installation load switch
- Quickly installs on left or right side of stairway
- Rack and pinion drive with concealed gear rack
- Low-to-the-floor footrest
- 3" to 6" horizontal extension at top landing
- 25 to 45 degree angle applications
- ETL listed ASME A17.5 / B44 and A18.1
- 350 lbs capacity
- Overspeed governor



85 degree
swivel seat
at top landing

Staying Home

OUTDOOR LEGACY STAIRLIFT

- Lowest priced outdoor unit
- Plated rack & pinion gear
- Anodized aluminum rail
- Protective coating on circuit board
- Plated chassis frame
- Stainless steel hardware
- Heavy duty canvas cover
- Sealed switches & waterproof connections
- Outdoor grade powder coating

OUTDOOR WARRANTY

1 Year warranty on BATTERIES

2 Year warranty on PARTS


3 Year warranty on DRIVE TRAIN & CIRCUIT BOARD

OUTDOOR PACKAGE INCLUDES HEAVY DUTY COVER
Cargo Box is optional equipment.



Technical Information

Dimensions subject to change without notice.



CONTACT US DIRECTLY AT:
816-380-2427



OR CALL US TOLL FREE AT:
877-378-4275

Staying Home Corp. • P.O. Box 37 • Harrisonville, MO 64701
2501 Anaconda Rd, Harrisonville, MO 64701
info@stayinghome.com • www.stayinghome.com