



Distributed exclusively by:



Otolift

Modul-Air Smart

Content

Redefining the Future of Mobility	1
Fast Delivery Curved Rail	3
Refreshing Elegance	5
Unique Rail System	7
Reverse Drive Technology	9
Flexible Convenience	11
Intuitive and Easy to Use	13
Rail Options	17
Key Benefits	19
Make It Match Your Decor	21



Otolift Modul-Air Smart

Redefining the Future of Mobility

Introducing Otolift's most advanced stairlift: the Modul-Air Smart. This stairlift pairs advanced safety technology with effortless operation, all wrapped in exceptional comfort and modern, high-quality upholstery for a sleek, stylish ride.

Every home, like its resident, is unique - that's why we designed our lifts to blend into your home decor. With the slim rail and narrow chair fold, this lift is a smart addition to your independent living journey. Paired with the innovative technology like silent, steady travel and display monitor on the chair arm, you'll glide swiftly into the future of safe home accessibility.



Otolift Modul-Air Smart

Fast Delivery Curved Rail

As the home access industry booms, so does the product technology. The Otolift's modular rail means that all rail pieces, sizes, and shapes are stockable - this way, you can enjoy your whole home safely and confidently within 10 days!



Otolift Modul-Air Smart

Refreshing Elegance

We introduced the elegant, inconspicuous design, but what about stairway use when the lift is parked?

When the chair is at the top or bottom of the steps, fold the seat, arms, and footrest up for optimal access to your stairs and walkway.

When in use, the chair silently travels up and down the stairs, keeping those midnight snack trips to the kitchen a secret. Use the control on the chair's arm to glide up and down to your next destination.



Otolift Modul-Air Smart

Unique Rail System


The Otolift Modul-Air Smart stands out in the curved stairlift category with its innovative modular rail system. Using a series of precisely engineered rail sections, it adapts seamlessly to nearly any curved staircase—delivering a tailored fit without the long wait times of fully custom-fabricated modular rails.



World's thinnest single rail
With a diameter of only 2 3⁄8" (6 cm) the rail of the Modul-Air Smart is the thinnest single rail in the world.



Unique angled coupling
Each rail component fits together seamlessly, creating a smooth, easy glide from start to end.



Tested for long-term use
The stairlift is thoroughly tested for a 10-year service life under constant heavy load. This ensures you can rely on long-lasting and dependable performance.



Otolift Modul-Air Smart
ø 2 3⁄8" (6 cm)
In full scale



Traditional stairlift
ø 3 1⁄8" (8 cm)
In full scale

Otolift Modul-Air Smart **Reverse Drive Technology**

Thanks to the optional power swivel upgrade, the chair can partially or even completely drive in reverse. Now you don't have to worry about squeezing past those tight turns and cramped spaces.



Otolift Modul-Air Smart

Flexible Convenience

The Otolift Modul-Air Smart is ideal for installation along the inside or outside of your stairs. Not only will you enjoy the powered folding footrest, but the swivel seat and smart dashboard in the chair's arm for stairlift diagnostics are great, intuitive additions for safety, comfort, and ease of use.

Wide, narrow or spiral staircases, the Otolift Modul-Air Smart easily adapts to the most demanding stairs and unobtrusively follows all curves. In addition, the rail takes up so little space that the stairs also can easily be used on foot.



Powered folding footrest

Powered footrest comes standard on all units. This means you don't have to bend down to fold the footrest up and down.



Optional powered swivel seat

The powered swivel upgrade allows you to control the swivel with the joystick controls. This also allows you to take advantage of the reverse drive technology on narrow staircases.



Smart dashboard

The smart dashboard in the armrest is a unique feature. It offers you a live view with the status of your stairlift. Simply check the battery level at a glance.



Otolift Modul-Air Smart

Intuitive and Easy to Use

The Otolift Modul-Air Smart is very easy to use. After installation, you will receive a calm, clear and detailed explanation from the installer.



Step 1

You bring the chair to you with the wireless call/send control.



Step 2

You sit down on the comfortable seat.



Step 3

The chair moves silently and smoothly.



Step 4

The chair continues its way up silently and stops automatically.



Step 5

The chair turns to the landing. Now you can take off your seat belt.



Step 6

You safely step off your stairlift.

Otolift Modul-Air Smart

Intuitive and Easy to Use



Smart battery system

Thanks to the smart battery system, the stairlift can also be used in the event of a power failure.



Effortless

The seat has ergonomic controls and is easy to use. The stairlift operates effortlessly, even with limited hand strength.



Simple

The new seat belt is easy to use and can be closed with one hand.



Safe

Your Otolift stops when it detects an obstacle. Automatic, gentle and safe.



Space-saving

The Otolift can be folded quickly. It saves space and the stairs remain easily accessible.



Quiet

The stairlift moves almost silently, producing less than 50 dBA, which is comparable to the sound level of a quiet dishwasher.

Adjust to situation

Rail Options

Find the perfect configuration for your staircase. With tailored start options and beautifully matched finishes, every choice ensures a smooth fit within your home. Whether installed on the inside or outside bend, the Otolift Modul-Air Smart adapts effortlessly to your living environment.

Start options

Explore a range of options for rail starts tailored to precisely meet each staircase.



Mint Start



Drop Start




Horizontal Start



Incline Start




Incline Start 90°




Incline Start 180°

Finish options


Discover a variety of options for rail finishes designed to perfectly match each unique situation.



Horizontal Finish



Horizontal Finish 90°



Horizontal Finish 180°

Inside or outside bend

Our stairlifts can be installed on either side of the staircase, optimizing space for safe passage by placing the Otolift Modul-Air Smart on the narrower side.



Inside Curve (left)



Outside Curve (right)

Otolift Modul-Air Smart

Key Benefits

- 1 Minimize handrail modifications
- 2 Customizable through an innovative modular rail system
- 3 The world's thinnest single rail
- 4 Smart dashboard in armrest
- 5 Ergonomic seat
- 6 Powered swivel seat option available
- 7 Powered folding footrest is standard
- 8 Rail supports mounted directly to your steps



Otolift Modul-Air Smart

Make It Match Your Decor

The seat of the Otolift Modul-Air Smart is made of high-quality synthetic leather and is available in two elegant colors: cream white (standard) and brown (optional).

The rail is supplied as standard in a sleek white finish that blends seamlessly into any interior.



The colors may vary depending on the print

Synthetic leather

-  Cream white (standard)
-  Brown (optional)







1500 Independence Blvd., Ste 220 | Sarasota, FL 34234
877-331-5438
otolift@harmar.com

[Harmar.com](https://www.harmar.com)